




The Naída Link CROS solution allows recipients to hear their best in quiet and in noise by providing sound awareness and ease of listening from their non-implanted side. The Naída Link CROS is designed for unilateral AB Naída CI recipients with an unaidable hearing loss in the non-implanted ear. The Naída Link CROS is worn on the contralateral ear and transmits the microphone input signal to the Naída CI processor via automatic, wireless, ear-to-ear network (HIBAN).



## Compatibility

Naída Link CROS is compatible with the Naída CI Q90 and Q70 sound processors.

## Components of the Naída Link CROS

1. Phonak CROS Demo Hook (not for daily usage)
2. Microphone input with microphone protector
3. Push button for mute functionality
4. Battery compartment (functions as ON/OFF switch)  
PowerOne size 13 batteries are required.



## Earpiece Options



### CROS SlimTube

- Recommended for everyday use
- Available in left/right and different lengths (size 0, 1, 2, 3)
- Domes: small/medium/large



### Customized CROS Tip

- Requires an earmold impression and is purchased from Phonak

## Setting up the Naída Link CROS Solution

- No fitting required to Naída Link CROS device
- Initialize the Naída CI processor in SoundWave™ 3.1 or later for use with the CROS device



### CROS Initialization

Select Right + CROS Left or Left + CROS Right




CROS icon defaults "ON" after CROS Initialization



## RogerReady


Automatically manage input from Naída CI microphones and Roger 17 receiver



- RogerReady is available in any Program Slot.
- RogerReady enabled by default for programs that utilize 50/50 or 30/70 audio mixing.
- When the Roger 17 receives the Roger input signal, the audio mixing will automatically change from Mic Only to 50/50 or 30/70 (as specified)

## Automated Progressive Programs

Automatically applies global M-level adjustment in steps of 5, 10, or 15 CU



- Simply set up the first Program Slot
- Then select 'Create Progressive Programs' from the slot drop down menu
- Choose the preferred increment
- New programs are automatically generated, saved in list of programs, and written to available Program Slots

## Custom Program Defaults

Select preferred program default settings based on the recipient's age



- Age-based default settings can be selected by the clinician under Preferences
- Create default settings for 5 different age groups
- Defaults are automatically applied to new programs, but do not alter copies of existing programs.