

Guide for Clinicians

Adult Cochlear Implant Home-Based Hearing Therapy



Cochlear®

Contents

Part 1: Purpose of the Home-Based Hearing Therapy Manual	4
Part 2: How to Use the Home-Based Hearing Therapy Manual	6
Part 3: Screening tests	10
1: Screening tests Levels A–D	10
2: Phoneme recognition tests	13

Appendix 1: Screening test forms A–D

Screening test forms	16
----------------------------	----

Appendix 2: Vowel and consonant confusion matrix forms

Matrix forms	22
--------------------	----

Authors	24
---------------	----

Acknowledgments	26
-----------------------	----

Part 1: Purpose of the Home-Based Hearing Therapy Manual

Hearing Therapy is an important component of adult post-implant rehabilitation. It can be useful to train listening skills and communication strategies, empower the recipient in their listening progress, and encourage good speaking habits in the family and with close associates. Also, areas of ongoing difficulty can be used as an indicator for further improvements of cochlear implant programming.

As clinicians find themselves managing ever-growing cohorts of recipients, it may be difficult to spend sufficient time with recipients with recently activated devices, training on auditory skills and building their confidence in the new auditory stimulus. Clinician-guided home practice consolidates and supplements the work done in the clinic and is an important component of the post-implant process in its own right. It enables the recipient to listen within their own environment, which includes such effects as reverberation and environmental distractions, rather than the controlled environmental conditions of the clinic room. In addition, the recipient has the opportunity to become more proficient at listening to their most frequent communication partners in their daily lives.

The Adult Cochlear Implant Home-Based Hearing Therapy Manual has been designed to provide hands-on materials that are flexible and adaptable enough to cover a range of auditory abilities. Many exercises have a component for recipients and families to create exercises on topics/vocabulary of relevance and interest to them.

The aim of the manual is to supplement the work of the clinician in the clinic, rather than replace it. The manual provides a Hearing Therapy Exercise Plan, which can be completed by the clinician and the recipient. The screening tests (see Appendix 1 and 2) will help establish an appropriate training level for each recipient.

Part 2: How to Use the Home-Based Hearing Therapy Manual

Where to start in the manual

The manual for recipients is structured into two sections: Section A “Getting started” and Section B “Making it more challenging”. It is structured according to the auditory development hierarchy.

- Manual for recipients with postlingual hearing loss (Postlingual)
 - **Section A:** Getting started. Sound and speech identification (closed set) and easier speech recognition (open set) exercises.
 - **Section B:** Making it more challenging. Open set speech recognition and comprehension exercises in quiet and background noise, and telephone training.

It is important that training be conducted at an appropriate level for the individual recipients' auditory skills, where the exercises are challenging but not too difficult. Training is more effective if the recipient works at a level where some success is achieved (70–90%), but not 100% proficiency. This maintains the recipient's motivation but leaves room for improvement. The screening tests provide an indication of the level at which to start. This is a guide only; the recipients' success with the listening tasks attempted in the clinic is your best indication. We recommend that all recipients complete the sound awareness environmental sounds checklist to start the process of learning to filter out unwanted sounds.

Tips to maximize success

During the clinic session, make a written plan for the exercises to be completed at home by completing the “Hearing Therapy Exercise Plan”.

- Suggest the appropriate training environment, for example, the quietest place in the house to start, and then as they need more challenging environments, the porch, family room with television on in the background, or in a room the recipient has reported as especially challenging.
- Ask the recipient to complete the recipient notes section when doing the planned exercises at home.
- Encourage the recipient and partner to work at a distance of about three feet initially, and find a quiet part of the home, to maximize early success.
- Suggest that the training is done at the same time every day, about five days each week for 20–30 minutes at a time. A regular time slot increases the ease with which they can incorporate hearing therapy into everyday life.

How to demonstrate the manual to your recipient and their speaking partner

- Using the prompts at the start of the exercise, explain the purpose and method of each task.
- Demonstrate the exercise using the first one or two items in the exercise. You may need to discuss the volume of speech you are using, and the speed of speech you feel is appropriate for the recipient.
- Demonstrate how to rehearse the target items while pointing to them.
- Show the recipient and speaking partner how to remove visual cues without holding a hand close to the mouth (recipient looks down, speaker sits beside the recipient on the implanted side).
- Demonstrate appropriate feedback; discourage the use of **“no”**, and demonstrate alternatives such as **“listen again”**, **“nearly right, listen to the end of the sentence again”**, or **“it sounds a bit like that word but it’s a different word”**.
- Discuss how to provide training and further cues (for example, some visual cues, word emphasis) when the recipient experiences difficulty.
- Ask the speaking partner to do the next item in the exercise. Check that the speaker is using an appropriate voice level/speed of delivery/feedback.
- Show the recipient how they can create more material—look for blank tables at the end of the exercise where the recipient and their family can create a similar task with a topic/vocabulary that is relevant to the recipient.
- Suggest other environments or particular family members for speaking partners.
- For some of the later exercises the recipient will work from their own sheets, which are available in the appendixes of the manual.

Progressing through the exercises

Explain to the speaking partner how to extend an exercise that is becoming too easy, using the following methods:

- Use longer words.
- Include less familiar vocabulary.
- Increase the distance (for example, by sitting across a table) or sit on the non-implanted side.
- Offer less repeats, encourage more guessing.
- Use less predictable vocabulary.
- Change the environment to incorporate progressively more noise.
- Slowly increase the rate of speech.

Part 3: Screening tests

It is important to use appropriate screening tests to assess the recipient's ability to identify and recognize speech. Firstly, the results from the screening tests may be used as a guide in setting an appropriate training level for the recipient. Secondly, repeated assessment using such screening tests before and after the training also provides some indications of whether recipients have benefited from the training or whether the training difficulty should be adjusted.

This manual provides four screening tests (Levels A-D). The screening forms are shown in Appendix 1 and can be photocopied. Screening tests Levels A and B are adapted from Mecklenburg, Dowell & Jenison (1982), Cochlear Adult Rehabilitation Manual.

1: Screening tests Levels A–D

Screening test Level A

- The recipient can look at the test sheet while you give the instructions.
- Instructions to recipient: **"I will say one word from each line, please say or point to the word I said."**
- Calculate the correct score as a percentage.
- Make a plan:
 - If score $\geq 80\%$: Give Screening test Level B.
 - If score $<80\%$: Have them start with Section A, Module 1 of the recipient's home-based hearing therapy manual.

Screening test Level B

- The recipient should look at the test sheet while you give the instructions.
- Instructions to recipient: **"I will say one of the four sentences in each set. Please point to or tell me which sentence I said."**
- Calculate the correct score as a percentage.
- Make a plan:
 - If score $\geq 80\%$: Give Screening test Level C.
 - If score $< 80\%$: Start at Section A, Module 4, of the recipient's home-based hearing therapy manual.

Screening test Level C

- The recipient must not see your copy of the test material. The recipient will have a copy with the first half of the sentence only.
- Instructions to recipient: **"I will say the first half of a sentence with you lip reading. I will then say the whole sentence without you lip reading. Please try to repeat the whole sentence. If you are not sure, please make the best guess that you can."**
- Score one point for each underlined word correctly repeated.
- Calculate the correct score as a percentage.
- Make a plan:
 - If score $\geq 80\%$: Give Screening test Level D.
 - If score $< 80\%$:
 - Use the phoneme confusion matrices (Appendix 2) to check vowel and consonant identification and specify training as required.

Screening test Level D

- The recipient does not see the test material as this is completely open set.
- Instructions to recipient: **"I will say some sentences. Please try to repeat as much of the sentence as you can, even if it's only one word. If you are not sure, make a guess. Guessing always counts."** Do not give repeats.
- Use the precise scoring method: Each word must be repeated exactly (for example, plurals are scored as incorrect—in sentence one, the response "dogs" instead of "dog" is scored as incorrect).
- Calculate the correct score as a percentage.
- Make a plan:
 - If score $\geq 80\%$:
 - Use the phoneme confusion matrices (Appendix 2) to check any residual vowel or consonant confusions, and give training as required.
 - Use the home-based hearing therapy manual: Start at Section B, Part 1, Module 8 and work through the manual from there (also include some of the earlier exercises in challenging conditions to enhance hearing in noise or over the phone).
 - If score $< 80\%$:
 - Use the phoneme confusion matrices (Appendix 2) to check any residual vowel or consonant confusions, and give training as required.
 - Use the home-based hearing therapy manual, Section A, Module 3 and work through the manual from there.

2: Phoneme recognition tests

Screening test Level A

It is important to reinforce the ability of recipients to identify phonemes. Vowel and consonant recognition is assessed to determine areas of difficulty then exercises are given based on the results. This guide includes phoneme recognition tests and includes phoneme identification exercises.

- The vowel recognition test (Appendix 2) requires the recipient to identify vowels from a set of 11 vowels (non-diphthongs), presented in an /h/–vowel–/d/ context (for example, "heed", "hid", "head" etc.). The consonant recognition test (Appendix 2) requires the recipient to identify consonants from a set of 12 consonants in an /a/–consonant–/a/ context (for example, "aba", "apa", "ama" etc.).
- The tables at the end of this section provide a list of the vowel and consonant exercises included in the manual.

Test procedures:

- Present the vowel or consonant recognition test by audition alone, randomly selecting stimulus items from the vowel matrix (for example, "heed", "hid", "head" etc.) or consonant matrix (for example, "aba", "apa", "ama" etc.) in Appendix 2. Do not give repeats.
- Ask the recipient to repeat the stimuli.
- Record the results in the matrix.
- Analyze the matrix for the errors/pattern of errors and determine which phonemes/phoneme features require specific hearing therapy.
- Make a plan: Based on the errors in phoneme recognition, select appropriate exercises from the home-based hearing therapy manual with reference to the Phoneme identification exercises on p13.
 - Example: You perform consonant recognition testing with a recipient with a postlingual hearing loss. Analysis of the matrix shows accurate identification of voicing and manner of articulation cues but some confusion, for example, among aba/ada/aga and among asa/asha/afa (that is, errors in identifying the place of articulation). In the home-based hearing therapy manual, concentrate on Module 2, Exercises 10–13.

Phoneme identification exercises

Home-Based Hearing Therapy Manual (Module 2)

Vowel identification – Level 1

- **Exercise 1:** Vowel duration

Vowel identification – Level 2

- **Exercise 2:** Vowel frequency cues (same duration)
- **Exercise 3:** Vowel frequency cues (same duration, differ F2 only)

Consonant identification – Level 1

- **Exercise 4:** Voicing

Consonant identification – Level 2

- **Exercise 5:** Manner of articulation – stops vs nasals (voiceless)
- **Exercise 6:** Manner of articulation – stops vs nasals (voiced)
- **Exercise 7:** Manner of articulation – fricatives vs nasals

Consonant identification – Level 3

- **Exercise 8:** Manner of articulation – lateral vs nasal
- **Exercise 9:** Manner of articulation – affricate vs fricative

Consonant identification – Level 4

- **Exercise 10:** Place of articulation – b, d, g
- **Exercise 11:** Place of articulation – p, t, k
- **Exercise 12:** Place of articulation – f, s, sh
- **Exercise 13:** Place of articulation – m, n

Appendix I: Screening test forms A–D

Screening Level A*

1	park	baseball	impossible
2	hot	ice cream	entertainment
3	boy	shipwreck	understanding
4	cow	rainbow	helicopter
5	duck	mushroom	watermelon
6	meat	cowboy	escalator
7	cold	footstep	demonstration
8	calm	workshop	television
9	cash	sandwich	qualification
10	tea	hotdog	difficulty
SCORE (keywords correct)			-----/10
			----- %

*Adapted from Mecklenburg, Dowell & Jenison (1982). Cochlear Adult Rehabilitation Manual.

Screening Level B*

1	The goldfish is in the bowl.	6	Come in the front door.
	There is no time.		There could be a mistake.
	Sunrise will be early.		Dogs bark too much.
	The girl has long hair.		We eat lunch at one o'clock.
	She has gone now.		I can see you.
2	Coffee has caffeine.	7	There's a full moon tonight.
	There is no sunshine today.		Watch the football game.
	The last flood was in March.		The farm is in the country.
	The teacher spoke too quietly.		You can travel by train.
	Please wait at home.		The red pen is on the desk.
3	He bought a new car.	8	She is ten years old.
	Drink orange juice for a cold.		Please serve the tea.
	The children went to bed.		It's a long distance call.
	We swim in the lake.		We need wood for the fireplace.
	Drive the car home.		The suitcase is brown.
4	Bring it here.	9	She has a fever.
	I'll be there in a minute.		The winter months are cold.
	Practice makes perfect.		Take your vitamins.
	A whale is a mammal.		Knock on the door.
	The dog is brown.		Eat fresh fruit.
SCORE (keywords correct)			-----/10
			----- %

*Adapted from Mecklenburg, Dowell & Jenison (1982). Cochlear Adult Rehabilitation Manual.

Screening Level C: Clinician copy

		score
1	The children played with... a <u>bat</u> and <u>ball</u> .	/2
2	Go to the bank to... <u>open</u> an <u>account</u> .	/2
3	He went shopping to buy... <u>milk</u> and <u>bread</u> .	/2
4	I'm looking for a house... with <u>three</u> bedrooms.	/2
5	She spent her vacation... <u>swimming</u> at the <u>beach</u> .	/2
6	He traveled to work... by <u>bus</u> and <u>train</u> .	/2
7	At the farm they saw... some <u>cows</u> and <u>horses</u> .	/2
8	Please wake me up... at <u>five o'clock</u> .	/2
9	She cleaned the house with... a <u>mop</u> and a <u>broom</u> .	/2
10	The weather will be... <u>warm</u> and <u>sunny</u> .	/2
	TOTAL (keywords correct)	-----/20
		----- %

Screening Level C: Recipient copy

1	The children played with...
2	Go to the bank to...
3	He went shopping to buy...
4	I'm looking for a house...
5	She spent her vacation...
6	He traveled to work...
7	At the farm they saw...
8	Please wake me up...
9	She cleaned the house with...
10	The weather will be...

Screening Level D*

		score
1	The <u>dog</u> barked at the <u>cat</u> .	/3
2	<u>Kick</u> the <u>ball</u> .	/2
3	<u>Wash</u> the <u>shirt</u> .	/2
4	She <u>added</u> the <u>sugar</u> .	/2
5	He <u>ordered</u> <u>apple pie</u> .	/3
6	She <u>sewed</u> the <u>coat</u> .	/2
7	The <u>play</u> is <u>starting soon</u> .	/3
8	There is a <u>horse</u> in the <u>paddock</u> .	/2
9	The <u>sun</u> <u>rose</u> in the <u>sky</u> .	/3
10	The <u>clothes</u> <u>dried outside</u> .	/3
	TOTAL (keywords correct)	-----/25
		----- %

*Adapted from IEEE sentences. Institute of Electrical and Electronic Engineers (1969).

Appendix 2: Vowel and consonant confusion matrix forms

Vowel matrix

Stimulus												Response											
	heed	hid	head	had	hard	hud	hod	hawed	hood	who'd	heard												
heed																							
hid																							
head																							
had																							
hard																							
hud																							
hod																							
hawed																							
hood																							
who'd																							
heard																							

Consonant matrix

Stimulus												Response											
	aba	apa	ama	ava	afa	ada	ata	ana	aza	asa	aga	aka											
aba																							
apa																							
ama																							
ava																							
afa																							
ada																							
ata																							
ana																							
aza																							
asa																							
aga																							
aka																							

Authors

Belinda Henry

PhD, DipAud, BSc (Hons), MAudSA (CCP)

Audiologist – Clinical Leader

Hearing Implant Program
Lady Cilento Children's Hospital
Queensland, Australia

Audiologist

Mater Cochlear Implant Clinic
Mater Health Services
Queensland, Australia

Academic Title Senior Lecturer

School of Health and Rehabilitation Sciences
University of Queensland
Queensland, Australia

Karen Pedley

MSc (Aud), BSc (Hons), MAudSA (CCP)

Audiologist

Manager, Attune Hearing Implant Centre
Queensland, Australia

Academic Title Senior Lecturer

School of Health and Rehabilitation Sciences
University of Queensland
Queensland, Australia

Qian-Jie Fu

PhD

Professor

Department of Head and Neck Surgery
David Geffen School of Medicine
University of California, Los Angeles
California, USA

Scientist III

House Research Institute
2100 West Third Street, Los Angeles
California, USA

Acknowledgments

Some of the exercises in the Adult Cochlear Implant Home-Based Hearing Therapy Manual are from the Cochlear Adult Auditory Rehabilitation Guide (Pedley, Lind & Hunt, 2005) and were originally included in the Cochlear Adult Rehabilitation Manual (Mecklenburg, Dowell & Jenison, 1982).

The authors wish to thank the cochlear implant recipients at the Mater Cochlear Implant Clinic, Queensland, Australia, and the Attune Hearing Implant Centre in Queensland, Australia, who trialled and provided valuable feedback on an earlier version of the manual, as well as Mater Health Services audiologists Barb Plath and Anika Batros, for providing valuable feedback on an earlier version of the manual and Guide for Clinicians.

Hear now. And always

As the global leader in implantable hearing solutions, Cochlear is dedicated to helping people with moderate to profound hearing loss experience a life full of hearing. We have provided more than 550,000 implantable devices, helping people of all ages to hear and connect with life's opportunities.

We aim to give people the best lifelong hearing experience and access to innovative future technologies. We have the industry's best clinical, research and support networks.

That's why more people choose Cochlear than any other hearing implant company.

©Cochlear Limited 2019. All rights reserved. Hear now. And always and other trademarks and registered trademarks are the property of Cochlear Limited or Cochlear Bone Anchored Solutions AB. The names of actual companies and products mentioned herein may be the trademarks of their respective owners.

www.Cochlear.com/US

Follow us on    

Cochlear Americas

13059 East Peakview Avenue
Centennial, CO 80111 USA
Telephone: 1 303 790 9010
Support: 1 800 483 3123

⋮ **Cochlear Canada Inc.**
2500-120 Adelaide Street West
Toronto, ON M5H 1T1 Canada
Support: 1 800 483 3123



FUN3569 ISS1 APR19

