



Building on the Success of Genesis AI

The Evolution of the Starkey
Genesis AI Product Family

STARKEY EDUCATION & TRAINING

GENESIS^{AI}

Webinar FAQs

How to Ask a Question

Click Q & A, type your question, click Send.

Where to get the Handout

A PDF handout for this presentation may be accessed through your Pending Courses on Audiology Online or by clicking on the link in Chat within the classroom.

How to Earn CEUs

CEUs are available with a paid Audiology Online Unlimited CEU Access membership. After completing the course, go to your Pending Courses on Audiology Online and pass the multiple-choice exam.

Assistance or Tech Support

- Call Audiology Online at 800-753-2160
- customerexperience@continued.com
- Use the Q & A window

Learning Outcomes

1

List at least three available hearing aid styles for the Genesis AI product family.

2

Identify the StarLink chargers compatible with Genesis AI rechargeable products.

3

Identify at least two ways Starkey's Pro8 HydraShield™ provides enhanced ingress protection to maximize durability and reliability.



Hear better. Live better.


Starkey Neuro Processor

INDUSTRY'S

MOST

ADVANCED

PROCESSOR TECHNOLOGY



Neuro Sound Technology

The **biggest leap** forward ever in hearing technology

Adapts more quickly to deliver crystal-clear,
transparent sound quality to patients



80
million
hearing adjustments
every hour


Processor + Battery Technology




Increased
Battery Capacity

Genesis AI Hearing Aids

With an expanded lineup of styles, Genesis AI delivers maximum comfort, quality, and discretion.



Going
Beyond
IP68



ALL
NEW

GENESIS^{AI}



Pro8 HydraShield™

8 points of protection



Recessed and Kevlar-reinforced receiver cable withstands daily handling.

Full silicone seal of case parts and charge contacts repels moisture.

Completely enclosed Li-Ion battery module creates another layer of protection from moisture ingress.


More robust circuit sealing and processes enhance durability.



Internal microphone mesh barrier provides final layer of defense.

Redesigned acoustic ports impede the entrance of foreign materials.

Internal and external protective coating surpasses its IP68 rating.

Continuous overlapping enclosure design provides a tighter seal.







Engineered from
the inside out to be
WATERPROOF



Humidity
Test



Salt Mist
Test



Rain Spray
Test



Targeted Sweat
Drip Test



Extended
Reliability Testing



Salt Brine
Immersion Test

IP68+

Styles




Standard

RIC RT

mRIC R

RIC 312



Custom

Wireless ITC|HS|ITE R

Wireless CIC

IIC|CIC NW

GENESIS^{AI} Design

Genesis AI introduces an
all-new RIC design crafted
from meticulous research to
ergonomically fit the ear



RIC RT

24 | 20 | 16 | 12



Up to
2X battery life
than our previous technology
on a single charge!

BATTERY LIFE

Up to
51 hours
on a single charge

Up to **29 hours**
(CROS)

Specifications

- Neuro Processor
- Dual Radio
- Sensor Technology
- Rechargeable
- Telecoil
- Rocker Switch
- IP 68
- CROS compatible
with Universal RIC RT CROS Transmitter

**Industry's longest lasting
rechargeable battery in a RIC**

mRIC R

24 | 20 | 16 | 12



20% smaller
with
over 65%
longer lasting battery!

BATTERY LIFE

Up to
41 hours
on a single charge

Specifications

- Neuro Processor
- 2.4 GHz
- Sensor Technology
- Push Button
- IP 68
- Rechargeable

**Industry's smallest
rechargeable RIC with
sensors!**

Specifications

- Neuro Processor and features
- Dual Radio - 2.4 GHz + NFMI
- Sensor Technology
- Rocker Switch
- IP 68
- CROS compatible
with Universal RIC 312 CROS Transmitter
- Size 312 zinc air battery
 - Up to 6 days of battery life
 - Up to 4 days (CROS)
- No Telecoil

ALL
NEW

RIC 312

24 | 20 | 16 | 12

BATTERY LIFE

Up to
6 days




Re-engineered Case

Revolutionary mic port openings that are integrated into the case and a barrier layer to further protect microphones from foreign material

Mic covers are no longer necessary!

Omnidirectional
Mic Opening

Directional
Mic Openings



All new everything. Including colors.

Modernized colors, updated to better match skin tones



Tech Black



Graphite
Gray



Silver



Beige



Caramel



Chestnut



White



SnapFit 2.0 Receiver*

featuring FlexForm™ Receiver Cable

New 10-pin connector

Reduction in cable diameter,
increased pliability, and Kevlar-
reinforcement for a discrete,
comfortable fit that's built to last

*Compatible with Genesis AI only

Updated Matrix Names

40 gain → L

50 gain → M

60 gain → P

70 gain → UP*

Absolute Power



Custom Cased

Available for All

*Custom Cased Required

Snap Fit 2.0

Receiver housing and cable lengths remain unchanged

- Existing Comfort Buds, earmolds, Hear Clear, and measuring tools are compatible

Black packaging



StarLink Chargers

Sleek, more modern design

Space for cables to maintain natural shape even with large earmolds

Universal USB-C port for charging

White LED lights

Magnetic cleaning brush for easy access

**Turbo charging available
on all StarLink chargers**



StarLink Chargers

StarLink Charger

Reimagined solution with
dexterity challenges in mind

StarLink Premium Mini

Small, highly portable design
with on-board battery

StarLink Charger
mRIC



StarLink Premium Mini
mRIC

StarLink Charger
RIC




StarLink Premium Mini
RIC

Charger Variants

RIC RT and mRIC R styles have unique chargers

- Regulatory marking under the lid helps identify charger
- Chargers are otherwise identical except for the charge port adaptors
- The “smaller” charge port adaptors are for the mRIC R



Laser Marking
Inside Lid



Charge Port
Adaptors

Charger LED Behavior

New, sleek white LED lights

Charger will display charge status of the hearing aid with the lowest battery level





LED Behavior

Hearing Aid Charging Times

Number of solid white lights	Hearing aid battery level	Hearing aid charging time (from full depletion)	
•	Up to 25%	7 mins	Turbo charge
• •	Up to 50%	1.5 hrs	Standard charge
• • •	Up to 75%	2.5 hrs	Standard charge
• • • •	100% charged	3.5 hrs	Standard charge



Ship Mode

Rechargeable RIC and custom hearing aids are delivered in “ship mode.”

Prior to first use:

StarLink Charger 2.0 must be plugged into a power source (e.g., a wall outlet)

Rechargeable RIC devices must be placed into the appropriate charger for 3 seconds and then removed

Rechargeable custom devices can be placed on the charger for 3 seconds and then removed or if there is a program button on the device, press the button to exit ship mode

Styles




Standard

RIC RT

mRIC R

RIC 312



Custom


Wireless ITC|HS|ITE R

Wireless CIC

IIC|CIC NW

2.4 GHz Custom R

24 | 20 | 16 | 12



BATTERY LIFE

Up to

42 hours

on a single charge

Specifications

- Neuro Processor
- 2.4 GHz
- Sensor Technology
- Telecoil (ITE only)
- IP 68
- Same industry-leading modular faceplate design as Evolv AI

2.4 GHz CIC

24 | 20 | 16 | 12



Industry's
smallest
wireless
custom

BATTERY LIFE

Up to

7 days

Specifications

- Neuro Processor and features
- 2.4 GHz with curved antenna (no NFMI)
- Sensor Technology
- Optional Push Button
- Tinnitus Technology
- IP 68
- Compatible with My Starkey App and StarLink Accessories
- 312 Zinc Air battery
- Up to 8 days battery life

CIC NW

24 | 20 | 16 | 12



BATTERY LIFE

Up to

7 days

Specifications

- Neuro Processor
- Optional Push Button
- IP 68
- Tinnitus Technology
- T2 Compatible
- Optional 312 battery
Up to 10 days

IIC NW

24 | 20 | 16 | 12



BATTERY LIFE

Up to

7 days

Specifications

- Neuro Processor
- Updated battery door design
- Tinnitus Technology
- IP 68
- 10 Battery
- T2 Compatible


**Industry's smallest
custom product**

Matrices



IIC NW

L- (35 gain), L, M



CIC

L, M, P, UP




**2.4 GHz Custom R
ITC**

L, M, P, UP



**2.4 GHz Custom R
ITE**

M, P, UP



GENESIS^{AI}

StarLink Accessories

Everything is better when you accessorize

Connect to favorite music, TV
shows, and phone calls with
superior sound quality

StarLink



*Legacy 2.4 GHz accessories are not compatible
with Genesis AI hearing aids*

StarLink

Remote Microphone +

Multipurpose, easy-to-use remote microphone

TV Streamer

Stream audio from a TV or other sources

Mini Remote Microphone

Miniature, single-purpose remote microphone

Table Microphone

Eight directional microphones and beam-forming technology enhances listening in challenging environments

Remote Control 2.0

Offers simple hearing aid control



ALL
NEW


StarLink Remote Control 2.0

Redesigned for maximum simplicity,
StarLink Remote Control 2.0 lets users easily
control volume, change programs, and
quickly access a favorite feature.



Remote Anatomy

StarLink Remote Control 2.0



Remote Anatomy

StarLink Remote Control 2.0

Field-replaceable cell battery (CR 2032)



Battery Door



User Controls

StarLink Remote Control 2.0

LED indicator

- White LED Illuminates for any button press

Volume Control

- Short press

Favorite Button

- May include: Home, Accessory Start/Stop
Edge Mode+*, Tinnitus On/Off, and
Jump to Program 4

**Edge Mode+ not available for 12 tier*



Accessory Pairing

StarLink Remote Control 2.0

- Power hearing aids off/on
- Press and hold the Favorite button for 5 seconds
- Indicator light becomes solid for 2 seconds for a successful pairing



Support Materials



Product Guide



Patient Brochure



Lifestyle Posters



Product and Wellness Posters



Replica Display


Starkeypro.com

White Papers

Quick Tips

Videos

Training Opportunities





GENESIS^{AI}

Anything But Ordinary.

