



BIMODAL FITTING REPORT



Pete Zahut

Birth Date: 03/06/1990

Ⓕ Cochlear Implant (CI)

Ⓓ Hearing Aid (HA)



Processor
Serial Number
Fitting Date

Naída CI Q90
1150495
09/01/2016



Phonak Naída Link RIC
or
Phonak Naída Link UP

PAIRING ID# 52722

Enter the pairing ID in Target™ to activate Binaural VoiceStream Technology™ between the cochlear implant and the hearing aid.

Then, configure program structure options manually.

Program Structure and Settings

HA Program	Program Name	Program Options	CI #
A	Automatic (Startup Program) <i>Startup beeps</i>	Default Settings	1
1	Acoustic Phone <i>2 beeps</i>	DuoPhone Preferred phone ear	2
2	StereoZoom <i>3 beeps</i>	Enabled Left	3



BIMODAL FITTING REPORT



Pete Zahut

Birth Date: 03/06/1990

Device Options

Settings

Push button	Program change (with QuickSync)
Volume button	Volume change (with QuickSync)
Volume behavior	Louder 10 steps, Softer 10 steps
Signals and alerts	Low battery - Enabled
	Volume control - Enabled
	Manual program change - Enabled

Accessories

ComPilot ComPilot is configured by the CI fitting software, SoundWave™. Do not connect or change the ComPilot configuration with the Target software.

Device Color (Please circle the CI color):

- | | | |
|---|---|---|
| <input type="radio"/> Sand Beige | <input checked="" type="radio"/> Petrol | <input checked="" type="radio"/> Caribbean Pirate |
| <input checked="" type="radio"/> Chestnut | <input type="radio"/> Alpine White | <input type="radio"/> Princess Pink |
| <input type="radio"/> Silver Gray | <input checked="" type="radio"/> Velvet Black | <input checked="" type="radio"/> Ruby |
| | <input type="radio"/> Other | |

! About this Report

This report provides the information that will allow your Phonak hearing aid to be configured so that the Binaural VoiceStream Technology between it and your Naída CI sound processor is activated. In order to enable this feature, the pairing ID listed must be entered in the Target software.

Note that entering the pairing ID will not automatically setup the program structure, the program options, or the device options. These must be setup by your hearing aid fitter in the hearing aid fitting software.

Programming the hearing aid will not affect the programming of the cochlear implant sound processor.

Reference the Naída Link Bimodal Desktop Fitting Guide for HA fitting details.