

# Phonak Naída™ Link CROS

## Product information



Phonak Naída Link CROS is a dedicated product for clients with an unaidable hearing loss in one ear. It is worn on the unaidable ear in combination with an Advanced Bionics Naída CI sound processor on the other ear. Naída Link CROS wirelessly transfers the microphone input signal to the Naída CI sound processor. The sound processor receives the signal which makes sounds from the unaidable side audible. Naída Link CROS has no acoustical output.

Naída Link CROS is only compatible to work with an Advanced Bionics Naída CI sound processor. It is a ready-to-wear solution that does not require fitting with Phonak Target or any other fitting software.

---

## Product description



Naída Link CROS

Battery size	13
Push button	•
Volume control	
Nano coating	•
Dimensions (L x W x D)	31.3 x 13.4 x 8.5 mm (1.23 x 0.53 x 0.33 in)
Weight	1.87 g (0,066 oz)

---

## Housing colors



P1

Sand Beige



P4

Chestnut



P6

Silver Gray



P8

Velvet Black



P9

Ruby



Q1

Petrol



T7

Alpine White

---

## Earpiece options

---

### CROS DemoHook

- Hook to hold CROS on the ear during fitting and demonstration
- For demo purposes only
- Included with any Phonak CROS device

---

### CROS SlimTube

- Standard SlimTube for CROS devices
- Available in left/right and four different lengths (0, 1, 2, 3)
- Domes need to be attached. Available as open, closed and power dome in three sizes (S, M, L)

---

### CROS Tip

- Custom made earpiece for CROS devices
  - Available in left/right and four different lengths (0, 1, 2, 3)
-

---

## Wireless Communication Portfolio

If a streaming signal is received by the AB Naída CI sound processor, the streaming signal from CROS is automatically put on hold.

The connection to Naída Link CROS will be automatically re-established when streaming from the accessory is stopped.

---

## Key features

### **Mute function**

Mutes the audio signal delivered from the Naída Link CROS. This allows for the convenient elimination of unwanted sound from the CROS side without the need to remove the device from the ear. This can help encourage all day use of the CROS device.

### **StereoZoom**

Zooms in on a single voice in a diffuse, noisy environment. StereoZoom is available only for use with the AB Naída CI Q90 sound processor.

### **Automatically matched microphone mode**

When the AB Naída CI sound processor switches (either manually or automatically) from omni microphone use to a beamformer, the Naída Link CROS will automatically match the microphone mode.

---

# Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we remain true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally.

[www.phonakpro.com](http://www.phonakpro.com)