

Manufacturer	Provider Guide	Patient Guide	Limitations	Phone Compatibility
GN ReSound App: ReSound Smart 3D	GNResound Assist Quick Guide	GNResound Remote Fine Tuning Quick Guide	No pairing of accessories <u>Asynchronous:</u> No feedback calibration; No receiver reconfiguration; 5 adjustment requests can be open and then they will not be able to submit any more until at least one is closed by the provider	GNResound App Phone Compatibility
Oticon App: Oticon RemoteCare	Oticon RemoteCare Quick Guide	Oticon Remote Care Instructions For Use	No changing acoustics; No feedback calibration; No Tinnitus Sound Generator adjustments; No in-situ audiometry	Oticon App Phone Compatibility

Manufacturer	Provider Guide	Patient Guide	Limitations	Phone Compatibility
Phonak App: MyPhonak	Phonak Remote Support Fitting Guide	Getting Started with Phonak Remote Support	No firmware updates; No re-setting trial period; No change in technology for trial HAs; Some limitations on MPO adjustments; No feedback calibration	Phonak App Phone Compatibility
Signia (<i>note: this information applies to Tele Care; not Tele Care Basic</i>) App: Signia App	Signia Tele Care Quick Guide	See the Provider Guide page labeled "Install the Signia App on your Patient's Smartphone"	No notch therapy testing	Signia App Phone Compatibility
Starkey App: Starkey Thrive Hearing Control	Starkey Hearing Care Anywhere Guide	Hearing Care Anywhere Requests	No adding or deleting programs; No firmware updates; No feedback calibration	Starkey App Phone Compatibility

Manufacturer	Provider Guide	Patient Guide	Limitations	Phone Compatibility
Unitron App: Unitron Remote Plus	Unitron Get to Know Remote Adjust	See the Provider Guide page labeled "Getting started with remote adjust"	No adding or deleting programs; No changing Bluetooth ears; Can't delete wireless hearing aids; No factory reset; No feedback calibration; No Tinnitus Sound Generator adjustments	Unitron App Phone Compatibility
Widex App: Widex Remote Care	Widex Remote Care Guide	Remote Visit How-To Guide for Widex	No programming limitations, although RemoteLink device is required	Widex App Phone Compatibility

From: Zitelli, L. (2021). 20Q: Remote hearing aid programming - getting started. *AudiologyOnline*, Article 27750. Available at www.audiologyonline.com

Information was current at time of publishing and is subject to change.